



Information Content of EXAFS (I)



Advanced Topics in
EXAFS Analysis

Bruce Ravel

Sometimes, we have *beautiful* data. This is the merge of 5 scans on a 50 nm film of GeSb on silica, at the Ge edge and measured in fluorescence at NSLS X23a2.



Here, I show a Fourier transform window of [3 : 13] and I suggest a fitting range of [1.7 : 4.7]. Applying the Nyquist criterion:

$$N_{idp} \approx \frac{2\Delta k \Delta R}{\pi} \approx 19$$

This gives us an upper bound of the information content of that portion of the EXAFS spectrum.

Introduction

Statistics

Information

Parameters

Other topics

Iffetit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples



Constraints

Restraints

Artemis

Information Content of EXAFS (II)

Sometimes, we have less-than-beautiful data. This is the merge of 42 scans on a solution containing 3 mM of Hg bound to a synthetic DNA complex, measured in fluorescence at APS 20BM.



Here, I show a Fourier transform window of [2 : 8.8] and I suggest a fitting range of [1 : 3]. Applying the Nyquist criterion:

$$N_{idp} \approx \frac{2\Delta k \Delta R}{\pi} \approx 8$$

This talk discusses strategies for dealing with severely limited information content.



Advanced Topics in
EXAFS Analysis

Bruce Ravel

Introduction

Statistics

Information

Parameters

Other topics

Iffffit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples

Constraints

Restraints

Artemis

What is This Nyquist Criterion?

Given that we apply Fourier analysis to $\chi(k)$, we can treat EXAFS as a signal processing problem. If

- The signal is ideally packed **and**
- The error in the fitting parameters is normally distributed **and**
- We understand and can enumerate all sources of error **and**
- We know the theoretical lineshape of our data **then**

$$N_{idp} \approx \frac{2\Delta k \Delta R}{\pi}$$

where, for EXAFS, Δk is the range of Fourier transform and ΔR is the range in R over which the fit is evaluated.

Unfortunately ...

None of those conditions really get met in EXAFS. N_{idp} is, at best, an upper bound of the actual information content of the EXAFS signal.



Advanced Topics in
EXAFS Analysis

Bruce Ravel

Introduction

Statistics

Information

Parameters

Other topics

lfeffit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples

Constraints

Restraints

Artemis

Statistical Parameters: Definitions



Advanced Topics in
EXAFS Analysis

Bruce Ravel

IFEFFIT uses a Levenberg-Marquardt non-linear least-squares minimization, a standard χ^2 fitting metric, and a simple definition of an R-factor:

$$\chi^2 = \frac{N_{idp}}{\epsilon N_{data}} \sum_{i=\min}^{\max} \left[\text{Re} (\chi_d(r_i) - \chi_t(r_i))^2 + \text{Im} (\chi_d(r_i) - \chi_t(r_i))^2 \right] \quad (1)$$

$$\chi_\nu^2 = \frac{\chi^2}{\nu} \quad (2)$$

$$\nu = N_{idp} - N_{var} \quad (3)$$

ϵ = measurement uncertainty

$$\mathcal{R} = \frac{\sum_{i=\min}^{\max} \left[\text{Re} (\chi_d(r_i) - \chi_t(r_i))^2 + \text{Im} (\chi_d(r_i) - \chi_t(r_i))^2 \right]}{\sum_{i=\min}^{\max} \left[\text{Re} (\chi_d(r_i))^2 + \text{Im} (\chi_d(r_i))^2 \right]} \quad (4)$$

Introduction

Statistics

Information

Parameters

Other topics

Iffeffit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples



Constraints

Restraints

Artemis

An Obviously Good Fit

Here is a fit to the first two shells of copper metal at 10 K



This is an unambiguously good fit:

\mathcal{R}	0.0025
N_{idp}	16
ν	12
S_0^2	0.95(3)
E_0	5.98(36) eV
a	3.6072(26) Å
Θ_D	505(16) K

Yet $\chi^2 = 32.03$!

What's goin' on here?

Why is χ^2 for an obviously good fit so much larger than 1?



Advanced Topics in
EXAFS Analysis

Bruce Ravel

Introduction

Statistics

Information

Parameters

Other topics

Iffit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples

Constraints

Restraints

Artemis

Statistical Parameters: Fit Evaluation



Advanced Topics in
EXAFS Analysis

Bruce Ravel

Introduction

Statistics
Information
Parameters
Other topics

Iffefit

MKW
Explain
Example

MFC
Explain
Example

MDS
Explain
Examples

Constraints

Restraints




Artemis

The determination of measurement uncertainty is, perhaps, a bit hokey in **IFEFFIT**. It is the average of the signal between 15 Å and 25 Å in the Fourier transform – a range that probably does not include much signal above the noise.

Is that signal between 15 Å and 25 Å in copper metal? Perhaps....

In any case, this method ignores the following:

- Approximations and errors in theory
- Sample inhomogeneity
- Detector nonlinearity
- Gremlins ;-)



Statistical Parameters: Interpretation I



Advanced Topics in
EXAFS Analysis

Bruce Ravel

OK then ... what is the implication of ϵ never being evaluated correctly by **IFEFFIT**?

- 1 χ^2_ν is always somewhere between big and enormous.
- 2 χ^2_ν is impossible to interpret for a *single* fit.
- 3 χ^2_ν **can** be used to compare different fits. A fit is improved if χ^2_ν is significantly smaller.
- 4 Error bars are taken from the diagonal of the covariance matrix. If χ^2_ν is way too big, the error bars will be way too small. The error bars reported by **IFEFFIT** have been scaled by $\sqrt{\chi^2_\nu}$.
- 5 Thus the error bars reported by **IFEFFIT** are of the “correct” size if we assume that the fit is a “good fit”.

Introduction

Statistics
Information
Parameters
Other topics

Iffeffit

MKW
Explain
Example

MFC
Explain
Example

MDS
Explain
Examples

Constraints

Restraints

Artemis

Statistical Parameters: Interpretation II



Advanced Topics in
EXAFS Analysis

Bruce Ravel

How do we know if a fit is “good”?

- The current fit is an improvement over the previous fit if χ^2_ν is sufficiently smaller.
- You should be suspicious of a fit for which N_{var} is close to N_{idp} , i.e. a fit for which ν is small.
- All variable parameters should have values that are physically defensible and error bars that make sense.
- The results should be consistent with other things you know about the sample.
- The R-factor should be small and the fit should closely overplot the data. (That was redundant. 😊)

Introduction

Statistics

Information

Parameters

Other topics

Ifeffit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples

Constraints

Restraints

Artemis

Interpreting Error Bars



Advanced Topics in
EXAFS Analysis

Bruce Ravel

Introduction

Statistics

Information

Parameters

Other topics

lfeffit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples

Constraints

Restraints

Artemis

The interpretation of an error bar depends on the meaning of the parameter.

- A fitted σ^2 value of, say, 0.00567 ± 0.00654 is troubling. That result means suggests that σ^2 is quite ill-determined for that path and not even positive definite. Yikes!
- On the other hand, a fitted E_0 value of, say, 0.12 ± 0.34 is just fine. E_0 can be positive or negative. A fitted value consistent with 0 suggests you chose E_0 wisely back in **ATHENA**.

Outside Knowledge

Because the information content of the XAS measurement is so limited, we are forced to incorporate knowledge from other measurements into our data analysis and its interpretation.

- Other XAS measurements — for instance, the “chemical transferability” of S_0^2
- Diffraction tells us structure, coordination number, bond lengths, etc.
- Things like NMR, UV/Vis, and IR can tell us about the ligation environment of the absorber
- Common sense:
 - $R_{NN} \not\approx 0.5 \text{ \AA}$, $R_{NN} \not\approx 4.0 \text{ \AA}$
 - $\sigma^2 \not\approx 0 \text{ \AA}^2$
- ... and anything else your (physical || chemical || biological || whatever) intuition tells you



Advanced Topics in
EXAFS Analysis

Bruce Ravel

Introduction

Statistics

Information

Parameters

Other topics

lfeffit

MKW

Explain

Example

MFC

Explain

Example

MDS

Explain

Examples

Constraints

Restraints

Artemis